

Together. All the way™

VitroTemp

For accurate, easy, and consistent temperature measurements of media during operations outside the incubator



Mastering complexity


Ensuring the right temperature during oocyte and embryo culture is critical for clinical success. However, minimising temperature deviations remains a technical challenge for many IVF labs, especially during handling outside the incubator. A dedicated solution for this purpose makes it easier to set an optimal and consistent temperature.

Temperature is critical for oocyte and embryo development

Ensuring that gametes and embryos are kept at appropriate temperatures during handling outside the incubator is a critical factor in the IVF lab workflow. Even small temperature deviations from 37°C can have a detrimental effect on oocyte and embryo development.

Temperatures above 37°C generally cause negative effects to oocytes and embryos. It is therefore very important to make sure the temperature in the media never exceeds 37°C. Embryos respond to heat through the expression of stress response genes. High temperatures can also result in conformation changes in molecules and structures, such as denaturation of proteins, which can have a detrimental effect on oocyte and embryo physiology and development¹.

Temperatures below 37°C can affect the stability of the meiotic spindle of oocytes, resulting in reduced fertilisation rates, delayed embryo development, and decreased clinical pregnancy rates^{2,3}.



“Controlling the medium temperature when dishes are handled outside the incubator is one of the most critical factors in the IVF lab workflow. This is of particular importance during lab procedures such as vitrification because the gametes and embryos are handled for longer periods outside the incubator. This solution is very helpful for that, as it perfectly reflects the temperature of the culture media surrounding the oocyte or embryo. We use it frequently, since it makes our work much easier. All our lab staff are very enthusiastic about it.”

Professor Thomas Freour, Chef de service – Biologie et Médecine de la Reproduction, gynécologie médicale, Responsable du centre AMP, CHU de Nantes, France

Meeting the challenge


Despite the use of heated surfaces, it is difficult to control media temperature in IVF dishes during manipulation outside the incubator. Measuring the temperature of heated surfaces alone is not sufficient, as media temperatures are affected by the interaction of both heated and non-heated surfaces and a range of other variables. Factors such as the temperature of the heated and non-heated surfaces as well as the ratios between these surfaces, evaporation, heat radiation, heat transfer, dish design, and dish setup all combine to affect media temperatures.

Temperature measurement of media in IVF dishes held on heated surfaces is further complicated by the actual process of temperature measurement. The properties of some temperature probes commonly used by IVF labs can affect their suitability for the various temperature measurements required.

A dedicated solution helps to overcome these issues, contributing to accurate, consistent, and reliable media temperatures.

“VitroTemp is the best temperature measurement equipment that we have introduced into our equipment maintenance program in the 30+ years that I have been in IVF. It provides us with a much greater level of confidence that the temperature we are recording is an accurate and true indication of what the embryo is being exposed to. It does not fluctuate as much as other probes, which has sped up our temperature checking processes in the morning. The probe is protected in the dish, so it is not prone to mechanical damage like so many others we have used, which means reduced replacement costs. I highly recommend that every IVF lab uses this setup. Our service contractors have also noted how much more accurate it is.”

*Dr Debbie Blake, PHD, Scientific / Quality Director
ReproMed, Auckland, New Zealand*



Optimise media temperatures in dishes with VitroTemp

VitroTemp™ makes it easier to maintain the optimal temperature of media in dishes based on accurate fast and simple measurements that replicate the actual conditions for procedures used in the IVF lab.


Take control of the dish temperature

VitroTemp is a specially designed temperature measurement solution that combines an advanced digital readout instrument with a custom made probe, a Pt100 sensor attached to the bottom of the Vitrolife square dishes. It's easy to use and provides more consistent temperature measurements than most thermometers commonly used in IVF labs today⁴.

With VitroTemp, it's much more straightforward to set the correct surface temperature for the Vitrolife dishes and media you will use when handling gametes and embryos outside the incubator. The measurement takes the exact dish and media setup into account, accurately replicating actual conditions for gametes and embryos.

How it works

Temperature measurement is based on electrical resistance in a small platinum sensor (Pt100) attached to the bottom of a Vitrolife square dish. Temperature is measured exclusively at this sensor, eliminating any thermal conductivity interference from the probe or the user. What's more, the reading is always accurate for the bottom of the dish, where gametes and embryos are positioned during handling and culture.



Built for the job

IVF labs put significant efforts into monitoring dish media temperatures. Unfortunately, not all temperature probes are suitable for this type of measurement due to their features and properties. All too frequently, labs use probes that are less than ideal for dish temperature measurement.

VitroTemp's purpose-designed probes help reduce the margin for error. Once you have set the optimal temperature for your chosen media and dish, it's easier to maintain the same temperature when you handle gametes and oocytes using the same size of Vitrolife dishes and the same media.

No more messy workarounds

VitroTemp eliminates the need to position and hold measuring equipment by hand, which can add complication and may interfere with measurements, causing errors and inconsistency. With VitroTemp, you just add medium (and oil if required), place the dish on the heated surface, and take your reading after allowing the medium to reach thermal equilibrium.

What's in the box?

VitroTemp comprises Vitrolife dishes in two sizes, each with a built-in Pt100 sensor, and a digital readout instrument.

Because the medium temperature depends on the dish type and setup, VitroTemp comes with two different types of Vitrolife dishes and probes to cover your needs for temperature calibration. The solution includes a 40 mm Culture dish and a 5 Well culture dish, each with a probe mounted in the bottom. An additional kit, including only the two probes and dishes, can be ordered separately.

Key benefits

- Accurate and sensitive
- Consistent optimal media temperature
- Simple to use
- Highly precise Pt-100 sensor
- Temperature measurement eliminates interference from thermal conductivity effects
- Temperature is always measured at the bottom of the dish, corresponding to embryo position



“Previously, we used a thin sensor that was dipped into medium drops or larger volumes in the dishes we used. It was critical that the sensor ended up correctly positioned in the medium. The measurements also had to be done without a lid. This gave an uncertain measurement. Now we have standardised the entire procedure using VitroTemp. With two types of dishes with integrated sensors, we can measure with the right setup for different workstations and work steps. We can also measure with or without a lid. Above all, VitroTemp has helped us standardise measurements and simplify the whole process.”

Maria Johansson, Biologist, Reproductive Medicine Sahlgrenska University Hospital, Gothenburg

Ordering information

Product	REF	Description
VitroTemp™	16210	Digital readout instrument and probe kit (5 Well + 40 mm)
VitroTemp™ probe kit	16211	Separate probe kit (5 Well + 40 mm dishes with built-in sensor)
VitroTemp™ calibration	16212	Calibration of probe, incl. shipping costs

Specifications: Testo 110 digital readout instrument

- Fast in-app configuration, graph history, second screen, and measurement data memory
- Quick and easy to use
- Highly precise

Protection class	IP40 (IP65 with optional TopSafe protective cover)
Power	3 x AA batteries
Battery life	100h

Important

VitroTemp dishes should never be in contact with gametes and embryos. This solution is strictly for temperature measurement and must not be used for holding or culturing gametes or embryos.

The dishes are in direct contact with the heated surface. This means that heated surface temperatures set for Vitrolife dishes will not be accurate for dishes without direct contact to the heated surface.

VitroTemp is specifically designed for temperature measurement of media in Vitrolife square dishes and should not be used with any other dishes.

The probes need to be calibrated on an annual basis. Find the list of accredited Vitrolife service partners here:



Square dishes specially designed for IVF

Vitrolife dishes are specifically designed to support and facilitate IVF procedures, making your work easier and more effective. All square dishes have an absolute flat bottom. When placed on a heated stage, all dishes receive the same bottom temperature. Once you have calibrated the temperature of your heated stage you can feel confident that all dishes will have the same temperature. Our dishes are manufactured with quality controlled and certified materials. The finished products are further tested to ensure a secure environment for gametes and embryos.



Micro-droplet culture dish

This unique innovation takes culture in micro-droplets to a higher level. The dish is equipped with 12 micro-wells optimised for 25-35 µL droplets.



5 Well culture dish

The 5 Well dish ensures embryo viability, as each well is surrounded by either air or media, resulting in a homogenous temperature.



Centre well dish

The Centre well dish is a multi-purpose dish for fertilisation, cryo-procedures, embryo culture, and embryo transfer. The oval centre well allows for easy instrument access from the sides. The small diameter of the bottom facilitates easy embryo location.



ICSI dish

The square format and low profile ICSI dish allows for easy instrument access.



Culture dish 40 mm

Like the other dishes intended for culture, this dish has tapered edges. The ramp makes for easier access to embryos as they are clearly visible when at the periphery of the well. The dish has a security ledge.

Guaranteed embryo and gamete safe

- Sterility assurance level (irradiation) 10^{-6}
- Less than 0.25 endotoxin units/device
- CE-marked for IVF
- MEA using multiple end-points, including 1-cell, expanded blastocyst within 96h $\geq 80\%$, and cell count

Product	REF	Description
Vitrolife square dishes		
Culture dish 40 mm	16001	For oocyte pick-up, fertilisation, and culture
Micro-droplet culture dish	16003	For embryo culture
5 Well Culture dish	16004	For denudation before ICSI, culture, and cryo procedures
Centre well dish	16005	For fertilisation, culture, cryo procedures, and transfer
ICSI dish	16006	For ICSI

Based on the regulatory status, not all products are available in all countries. Please contact your local sales representative for information on availability in your country.

Orders & customer support

Contact your local sales representative for prices and availability. Orders can be placed through our website at www.vitrolife.com. You can also contact us by email and phone:

Sweden office: Phone: +46 31 721 81 00 Email: order@vitrolife.com	Japan office: Phone: +81 03 6459 4437 Email: japan@vitrolife.com
US office: Phone: +1 866 848 7687 Email: order.us@vitrolife.com	China office: Phone: +86 10 6403 6613 Email: order.asia@vitrolife.com
Australia office: Phone: +61 1800 848 765 Email: order.australia@vitrolife.com	Germany office: Phone: +49 871 4306570 Email: germany@vitrolife.com
Denmark office: Phone: +45 7221 79 00 Email: order@vitrolife.com	

Technical support

Europe, Middle East & Africa: support.fertility@vitrolife.com
Americas: support.us.fertility@vitrolife.com
Asia: support.asia@vitrolife.com
Japan & Pacific: info@vitrolife.co.jp



Learn more about
VitroTemp and
Vitrolife square
dishes