EN: Indication for use

Media for immobilization and isolation of sperm prior to intracytoplasmic sperm injection, ICSI.

Product Description

ICSI™ is a viscous sperm handling solution containing polyvinylpyrrolidone (PVP) and recombinant human albumin (RHA).

ICSI™ is antibiotic free.

For use after equilibration at +20 ± 5°C and ambient atmosphere.

Storage instructions and stability

Store dark at +2 to +8°C.

ICSI™ is stable until the expiry date shown on the container labels and the LOT-specific Certificate of Analysis.

Media bottles should not be stored after opening. Discard excess media after completion of the procedure.

Directions for use

Prepare warmed and equilibrated dishes with droplets of handling-medium and a droplet of ICSI™ in the center.

Cover with OVOIL™.

Place the oocytes into the handling medium droplets, one oocyte per droplet.

Place a small volume (1–2 µL) of prepared sperm suspension into the center of the ICSI™-droplet. Incubate the dishes for a few minutes to allow the sperm to migrate to the outer perimeter of the droplet.

Prime the injection pipette with ICSI™.

Immobilize individual sperm by using the injection pipette to crush the membrane of the sperm tail. Aspirate the immobilized sperm.

Move the oocyte droplet into the field of view. Use a holding pipette to secure the oocyte for injection. Inject the sperm with the least amount of ICSI™. Place the injected oocytes in culture medium.

Specifications

Sterile filtered SAL 10³

Bacterial endotoxins (LAL assay) < 0.25 EU/mL

Human sperm motility recovery assay ≥40% of initial motility

Microscopic evaluation Passes test

Mouse Embryo Assay Not performed

LOT specific test results are available on the Certificate of Analysis provided with each delivery.

Precautions

Discard product if bottle integrity is compromised. Do not use ICSI™ if it appears cloudy.

ICSI™ contains polyvinylpyrrolidone and recombinant human albumin, but no antibiotics.

Re-use may result in microbiological contamination and/or property changes in the product.

To avoid contamination Vitrolife strongly recommends that media should be opened and used only with aseptic technique.

The risk of reproductive toxicity and developmental toxicity for IVF media, including Vitrolife's IVF media, have not been determined and are uncertain.

Not for injection.

Caution: Federal (US) law restricts this device to sale by or on the order of a physician.

Description of ISO Symbols

A Sterilized using aseptic processing techniques

Temperature limit

Do not re-use, discard after procedure.

Date of manufacture.

Use-by date.

Caution: Consult accompanying documents.

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