# OVOIL HEAVY™

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## CS: Indikace pro use

Oil for covering of medium during in vitro fertilisation and micro-manipulation procedures.

## Product Description

OVOIL HEAVY™ is a sterile, high viscosity oil. For use after equilibration at +37°C in CO2-incubator or ambient atmosphere, depending on the intended use.

OVOIL HEAVY is stable until the expiry date shown on the container labels and the LOT-specific Certificate of Analysis.

Bottles can be used for up to two weeks after first opening, use aseptic technique and minimize the time outside the refrigerator. Discard excess product no later than two weeks after first opening.

Using a sterile, non-toxic pipette, add medium to a non-toxic culture dish and cover entirely with OVOIL HEAVY.

Use aseptic technique and minimize the time outside the refrigerator. Discard excess product no later than two weeks after first opening.

Using a sterile, non-toxic pipette, add medium to a non-toxic culture dish and cover entirely with OVOIL HEAVY.

## Specifications

<p>Asceptically filtered</p> <p>Mouse Embryo Assay one-cell, expanded blastocyst within 96 hours [%]</p> <p>≥80</p> <p>Mouse Embryo Assay one-cell, cell number within 96 hours</p> <p>No statistical difference <sup>1)</sup></p> <p>Micro droplet Mouse Embryo Assay one-cell, expanded blastocyst within 96 hours [%]</p> <p>PASS</p> <p>Micro droplet Mouse Embryo Assay one-cell, cell number within 96 hours</p> <p>No statistical difference <sup>1)</sup></p> <p>Human Sperm Oil Assay</p> <p>PASS</p> <p>Bacterial endotoxins (LAL assay) [EU/mL]</p> <p>&lt; 0.25</p>	
---	--

<sup>1)</sup> Mean cell number of the test group is statistically compared to the control group. The mean of the test group must not be statistically different (P>0.05) to the control to constitute a Pass.

LOT specific test results are available on the Certificate of Analysis provided with each delivery.

## Precautions

Discard product if bottle integrity is compromised. Do not use OVOIL HEAVY if it appears cloudy.

To avoid contamination Vitrolife strongly recommends that bottle should be opened and used only with aseptic technique.

The risks of reproductive toxicity and developmental toxicity for IVF media, including Vitrolife’s IVF media, have not been determined and are uncertain.

Not for injection.

Caution: Federal (US) law restricts this device to sale by or on the order of a physician.

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