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The screenshot displays the 'Models' tab in the KIDScore software. On the left, there is a table listing models:

Active	Name	Type	Creator	Date
<input checked="" type="checkbox"/>	KIDScoreD5 v3.2	Imported	Vitrolife	2019-10-01
<input checked="" type="checkbox"/>	KIDScoreD3 v1.3	Imported	Vitrolife	2022-03-16

Below this table are 'Import' and 'Export' buttons. The 'Selected model' section shows a table of variables:

Variable	Min	Max
Information - NOT2PN	-	-
Information - 9NF	-	-
Information - t2	-	-
Information - t3	-	-
Information - t4	-	-
Information - t5	-	-
Information - t8	-	-
Information - cells66	-	-

The main area shows the detailed view for 'KIDScore D3 v1.3'. It includes fields for 'Model Name', 'Model Type' (Imported), and 'Creator' (Vitrolife). The 'Model Description' states: 'KIDScore D3 is defined by Vitrolife A/S based on the knowledge and experience extracted from our available KID data (please see the use manual for the EmbryoViewer software for a definition of KID data). The model focuses on which embryos to avoid rather than which embryos to select. It is thus a model which is based on avoidance criteria rather than selection criteria. The model will apply the avoidance criteria to the embryos and...'. Below this is a 'Custom Expressions' table:

Name	Expression
cells66	cells(66)

There are 'New', 'Edit', and 'Delete' buttons for custom expressions. The 'Model Definition' table is as follows:

Variable	Description	Min	Max	Classification
NOT2PN	Info			
9NF	Info			
t2	Info			
t3	Info			
t4	Info			
t5	Info			
t8	Info			
cells66	Info			

At the bottom of the detailed view are 'Save' and 'Clear' buttons. On the right side, there is a 'Model Provided By:' section with the Vitrolife logo and a license agreement section. The license agreement includes the text: 'The installation and use of this KIDScore D3 model (the "Model") shall be subject to the terms and conditions stated below. By clicking the "I agree" button or by installing or otherwise using the Model you have conclusively accepted to be bound by all of these terms and conditions. If you do not agree to all of the terms and conditions please do not install or use the Model. All rights in the Model belong to Vitrolife A/S ("Vitrolife"). Vitrolife grants you a non-exclusive, non-transferable, and non-sub-licensable license to install the Model and use it solely together with your EmbryoViewer and ES Server software with a valid license from Vitrolife. You are not granted any other rights or license with respect to the Model. Without limiting the foregoing you shall not copy, modify, decompile, reverse engineer, disassemble, or convert the Model or assign, transfer, sell, rent or lease the Model to any third party. Any actions, use, copying or distribution of the Model not authorized under these terms of use shall automatically terminate your rights hereunder. VITROLIFE DISCLAIMS ALL WARRANTIES EXPRESS OR IMPLIED, INCLUDING WITHOUT LIMITATION ANY IMPLIED WARRANTIES OF FITNESS FOR A PARTICULAR PURPOSE, MERCHANTABILITY OR NON-INFRINGEMENT OF THIRD-PARTY RIGHTS. VITROLIFE SHALL NOT UNDER ANY CIRCUMSTANCES BE LIABLE TO YOU OR ANY THIRD PARTY FOR ANY LOSS OF PROFITS, LOSS OF USE, INTERRUPTION OF BUSINESS, OR ANY INDIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES OF ANY KIND. This end-user license agreement shall be governed and interpreted solely in...'. At the bottom right, there is a 'KIDScore D3' logo with 'Software version 1', 'REF 15531', 'VERSION 1.3.0.27527', and 'CE MD' marks. The contact information 'UDI (01) 05712714675317 (0012) 1.3.0.27527' is also present.

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
The screenshot displays the KIDScore software interface. On the left, there is a large grayscale image of an embryo. Below the image is a progress bar with a red line indicating the current step. On the right, there is a table with columns: Anno. Tool, Value, Well, Dec., and Progress. The table lists various parameters and their values for well AA-1. Below the table, there is a 'Comments' field, 'Prev' and 'Next' buttons, and a 'Strategy' dropdown menu set to 'KIDScore D5 v3 Annotation'. At the bottom right, there are radio buttons for 'Step by' with 'Variable' selected and 'Well' unselected.

Anno. Tool	Value	Well	Dec.	Progress
PN	2	AA-1		100 %
L2	28.5 h	AA-2		0 %
L3	30.9 h	AA-3		0 %
L4	31.2 h	AA-4		0 %
L5	42.8 h	AA-5		0 %
IB	77.6 h	AA-6		0 %
ICM	A	AA-7		0 %
TE	B	AA-8		0 %
		AA-9		0 %
		AA-10		0 %
		AA-11		0 %
		AA-12		0 %
		AA-13		0 %
		AA-14		0 %
		AA-15		0 %
		AA-16		0 %

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Well	Dec.	Current score	NOT2PN	t2	t3	t4	t5	IB	ICM	TE	Last stage	Morph. grade	Last image	Saved score
AB-1		6.1	●	25.5	38.3	39.6	50.5	109.7	B	B		B		
AB-2		6.6	●	25.8	38.3	40.2	53.6	103.7	B	B		EB		
AB-3		8.4	●	27.7	40.5	41.3	53.9	103.9	B	A		EB		
AB-4		6.2	●	25.8	39.1	39.2	55.5	106.6	B	B		EB		

Current Model
 KIDScoreD5 v3
 Created 2018-11-01 by Vitrolife

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Well	Dec.	Current score	UNEVEN2	UNEVEN4	MN2	MN4	Frag-2	Coll. Count	Last stage	Morph. grade	Last image	Saved score
AB-1		NA	●	●	●	●	5.0	0.0		B		6.1
AB-2		NA	●	●	●	●	5.0	1.0		EB		6.6
AB-3		NA	●	●	●	●	20.0	0.0		EB		8.4
AB-4		NA	●	●	●	●	10.0	2.0		EB		6.2

Current Model
 Info example
 Created 2019-06-04 by ADMIN

Saved Model
 KIDScoreD5 v3
 Saved 2019-07-02 15:57:22 by ADMIN

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©2022 Vitrolife A/S. 保留所有權利。