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| Active | Name | Type | Creator | Date |
|-------------------------------------|------------|----------|---------|------------|
| <input checked="" type="checkbox"/> | Model name | Additive | ADMIN | 2014-10-28 |

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The screenshot displays the 'Models' tab in the KIDScore software. It features a table of models, a 'Selected model' section, and a detailed view of the 'KIDScore D3 v1.3' model. The detailed view includes the model name, description, creator (Vitrolife), and a table of variables with their descriptions and classifications. The license agreement is also visible at the bottom right.

| Active | Name | Type | Creator | Date |
|-------------------------------------|-----------------|----------|-----------|------------|
| <input checked="" type="checkbox"/> | KIDScoreD5 v3.2 | Imported | Vitrolife | 2019-10-01 |
| <input checked="" type="checkbox"/> | KIDScoreD3 v1.3 | Imported | Vitrolife | 2022-03-16 |

Selected model

| Variable | Min | Max |
|-----------------------|-----|-----|
| Information - NOT2PN | - | - |
| Information - IPNF | - | - |
| Information - t2 | - | - |
| Information - t3 | - | - |
| Information - t4 | - | - |
| Information - t5 | - | - |
| Information - t8 | - | - |
| Information - cells66 | - | - |

Model Definition

| Variable | Description | Min | Max | Classification |
|----------|-------------|-----|-----|----------------|
| NOT2PN | Info | | | |
| IPNF | Info | | | |
| t2 | Info | | | |
| t3 | Info | | | |
| t4 | Info | | | |
| t5 | Info | | | |
| t8 | Info | | | |
| cells66 | Info | | | |

Model Description

KIDScore D3 is defined by Vitrolife A/S based on the knowledge and experience extracted from our available KID data (please see the use manual for the EmbryoViewer software for a definition of KID data).

The model focuses on which embryos to avoid rather than which embryos to select. It is thus a model which is based on avoidance criteria rather than selection criteria. The model will avoid the avoidance criteria in the embryo and

Model Provided By:

Vitrolife

Vitrolife A/S End-User License Agreement for KIDScore D3 Model

The installation and use of this KIDScore D3 model (the "Model") shall be subject to the terms and conditions stated below. By clicking the "I agree" button or by installing or otherwise using the Model you have conclusively accepted to be bound by all of these terms and conditions. If you do not agree to all of the terms and conditions please do not install or use the Model.

All rights in the Model belong to Vitrolife A/S ("Vitrolife"). Vitrolife grants you a non-exclusive, non-transferable, and non-sub-licensable license to install the Model and use it solely together with your EmbryoViewer and ES Server software with a valid license from Vitrolife.

You are not granted any other rights or license with respect to the Model. Without limiting the foregoing you shall not copy, modify, decompile, reverse engineer, disassemble, or convert the Model or assign, transfer, sell, rent or lease the Model to any third party. Any actions, use, copying or distribution of the Model not authorized under these terms of use shall automatically terminate your rights hereunder.

VITROLIFE DISCLAIMS ALL WARRANTIES EXPRESS OR IMPLIED, INCLUDING WITHOUT LIMITATION ANY IMPLIED WARRANTIES OF FITNESS FOR A PARTICULAR PURPOSE, MERCHANTABILITY OR NONINFRINGEMENT OF THIRD PARTY RIGHTS. VITROLIFE SHALL NOT UNDER ANY CIRCUMSTANCES BE LIABLE TO YOU OR ANY THIRD PARTY FOR ANY LOSS OF PROFITS, LOSS OF USE, INTERRUPTION OF BUSINESS, OR ANY INDIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES OF ANY KIND.

This end-user license agreement shall be governed and interpreted solely in

KIDScore D3
Software version 1
[REF] 15531
VERSION 1.3.0.27527

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| Well | Dec. | Current score | NOT2PN | t2 | t3 | t4 | t5 | tB | ICM | TE | Last stage | Morph. grade | Last image | Saved score |
|------|------|---------------|--------|------|------|------|------|-------|-----|----|------------|--------------|------------|-------------|
| AB-1 | | 6.1 | ● | 25.5 | 38.3 | 39.6 | 50.5 | 109.7 | B | B | | B | | |
| AB-2 | | 6.6 | ● | 25.8 | 38.3 | 40.2 | 53.6 | 103.7 | B | B | | EB | | |
| AB-3 | | 8.4 | ● | 27.7 | 40.5 | 41.3 | 53.9 | 103.9 | B | A | | EB | | |
| AB-4 | | 6.2 | ● | 25.8 | 39.1 | 39.2 | 55.5 | 106.6 | B | B | | EB | | |

Current Model
 KIDScoreB5 v3
 Created 2018-11-01 by Vitrolife

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| Well | Dec. | Current score | UNEVEN2 | UNEVEN4 | MN2 | MN4 | Frag-2 | Coll. Count | Last stage | Morph. grade | Last image | Saved score |
|------|------|---------------|---------|---------|-----|-----|--------|-------------|------------|--------------|------------|-------------|
| AB-1 | | NA | ● | ● | ● | ● | 5.0 | 0.0 | | B | | 6.1 |
| AB-2 | | NA | ● | ● | ● | ● | 5.0 | 1.0 | | EB | | 6.6 |
| AB-3 | | NA | ● | ● | ● | ● | 20.0 | 0.0 | | EB | | 8.4 |
| AB-4 | | NA | ● | ● | ● | ● | 10.0 | 2.0 | | EB | | 6.2 |

Current Model
 Info example
 Created 2019-06-04 by ADMIN

Saved Model
 KIDScoreB5 v3
 Saved 2019-07-02 15:57:22 by ADMIN

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©2022 Vitrolife A/S. Visas tiesības aizsargātas.
VITROLIFE A/S, DĀNIJA