

Introducing EmbryoSlide+ ic8 dish. Expanding your options.

By popular demand, we now offer the possibility to culture embryos individually in culture wells in the EmbryoScope+ and EmbryoScope 8. Expanding options to fit your handling needs.

New



- **Uniquely designed culture dish for use in EmbryoScope+ and EmbryoScope 8 time-lapse instruments**
- **Designated area for barcode label for patient traceability**
- **8 individual wells for embryo culture and 4 rinsing wells**
- **Media and oil usage:**
20 µl media for culture well
25-30 µl media for rinsing well
1.6 ml oil overlay
- **Capacity of 8 embryos**

Orders & customer support

Ordering order@vitrolife.com
Global tech support support.embryoscope@vitrolife.com
+45 7023 0500 (24-hour hotline outside US)
888-879-9092 (US support)

EmbryoSlide+ ic8 dish has not received FDA 510k clearance.

Vitrolife A/S • Jens Juuls Vej 16 • 8260 Viby J • Denmark • Tel +45 7221 7900 • order@vitrolife.com • vitrolife.com

Vitrolife 
PART OF VITROLIFE GROUP