

Together. All the way™

Introducing EmbryoSlide+ ic8 dish. Expanding your options.

By popular demand, we now offer the possibility to culture embryos individually in the EmbryoScope+.

New



- Uniquely designed culture dish for use in EmbryoScope+
- Designated area for barcode label for patient traceability
- 8 individual wells for embryo culture and 4 rinsing wells
- Media and oil usage:
 - 20 µl media for culture well
 - 25-30 µl media for rinsing well
 - 1.6 ml oil overlay
- Capacity of 8 embryos

Orders & customer support

Ordering order@vitrolife.com
Global tech support support.embryoscope@vitrolife.com
+45 7023 0500 (24-hour hotline outside US)
888-879-9092 (US support)

Vitrolife 
PART OF VITROLIFE GROUP